

# Sirafan Speed RAPID DISINFECTANT

**Sirafan Speed** is a ready-to-use no rinse rapid surface disinfectant with a broad efficacy against bacteria and enveloped viruses. The synergistic action of its active ingredients provides quick acting disinfection on all food contact surfaces such as preparation tables, cutters and slicers. No rinse needed.

Bactericidal acc. EN 1276, Yeasticidal acc. EN 1650, Bactericidal and Yeasticidal acc. EN 13697, effective against Adenovirus + Murine Norovirus.

**Application area:** Food preparation surfaces, counters, buffets, kitchen equipment and utensils, hand contact surfaces like door handles, switches, etc.


- **Performance:** Leaves shiny surfaces and without smears.
- Efficiency: Fast disinfection action. Rinse free.
- Safety: Broad disinfection spectrum, prevents cross contamination and improves food safety in Retail and Food Service Industry.
























- 1 Check hygiene plan and product label for detailed use information.
- 2 Flammable product, do not spray on a flame or any incandescent object or hot surfaces.
- 3 Unplug all the electric equipment before use.
- 4 Wear suitable protective clothing when handling the product.
- 5 Spray on the surface to be disinfected.
- 6 Wipe with a clean and dry cloth.
- 7 Allow the product work for the indicated contact time:
  - EN 1276, 30 sec.
  - EN 1650, 30 sec.
  - EN 13697, 30 sec.
  - Adenovirus + Murine Norovirus - 10 min.

8 Leave to air dry. No rinse necessary.

## SPECIAL NOTES / INFORMATION

100 g of liquid contains: 35 g 2-propanol, 25 g 1-propanol **Use disinfectants safely. Always read the label and product information before use.** Dispose of packaging only if empty and closed.

For professional use only.

#### SAFETY INFORMATION

Avoid all possible sources of ignition. For further information and safety instructions please see product label and Safety Data Sheet.

#### STORAGE CONDITIONS

Store between 0 °C and 25 °C in original closed container only.

### PACK SIZE

6 x 750 ml bottles with trigger sprayer 2 x 5 L  $\,$ 





Specialist Ecolab Institutional × (+373) 79 101 029 www.toreco.md

